

A manual autoclave that sets the standard for steam sterilization

The Tuttnauer 2540M is a manual model with durable components that reduce running and maintenance costs. Tuttnauer's M Series Autoclaves are affordable options for clinics that do not want to compromise on quality, safety and reliability.

Features:

- Automatic shut off at the end of sterilization and dry cycles
- Electro polished chamber and door
- Double safety locking device prevents door from opening while chamber is pressurized
- Front-located drain valve for quick and easy draining
- Dual safety thermostat protects against overheating
- 2 year warranty



Specifications:

- Chamber dimensions: 10" (254 mm) x 19" (475 mm)
- Volume: 6 gallons | 23 Liters
- Overall dimensions: 20" (510 mm) W x 21.5" (545 mm) D x 14.4" (365 mm) H
- Standard cassette capacity: 3 full and 3 half
- Tray dimensions: 6.7" (170 mm) W x 16.3" (415 mm) x 0.8" (20 mm) H
- Number of trays: 4
- Cold cycle time: 30 minutes
- Hot cycle time: 14 minutes
- Voltage: 120V
- Frequency: 50/60 Hz
- Power: 1400W
- Current: 12A
- Shipping weight: 85 lbs. (39 kg)